

ZENO[®]

HOT SPOT

QUICK START GUIDE

HOW TO USE ZENO



STEP 1

Press the Power Button. The power indicator will light up and you will hear a tone to signal the treatment has started.



STEP 2

Place and hold tip on your blemish for the 2½ minute treatment cycle. You will hear a single tone every 30 seconds.

STEP 3

Remove ZENO from acne blemish when you hear a multi-tone sound, indicating the treatment is over.

Uses Remaining Indicator



60-80
uses
remain



40-59
uses
remain



20-39
uses
remain



0-19
uses
remain

TY4032-001